



# Comparison of Actual Performance in the Flow and Fraction of Inspired O<sub>2</sub> among Different High-Flow Nasal Cannula Devices: A Bench Study

**Authors:** Yuyan Zhou, Zhong Ni, Yuenan Ni, Binmiao Liang, and Zongan Liang

Department of Respiratory and Critical Care Medicine, West China Hospital, Sichuan University, No. 37 Guoxue Alley, Chengdu 610041, Sichuan, China

## Objectives


High-flow nasal cannula (HFNC) oxygen therapy has been recommended for use in patients with acute respiratory failure and many other clinical conditions. HFNC devices produced by different manufacturers may have varied performance. Whether there is a difference in these devices and the extent of the differences in performance remain unknown. This study sought to compare flow and FiO<sub>2</sub> at the cannula among different HFNC systems.

## Study methods

This was a bench study comparing four HFNC devices (AIRVO 2, TNI softFlow 50, HUMID-BH, and OH-70C) and one ventilator with a HFNC module (bellavista 1000). A gas flow analyzer was used to measure actual flow rate of gas delivered. An oxygen analyzer measured FiO<sub>2</sub>. High flow oxygen was provided by the hospital intensive care unit systems. Each device was measured at flow rates ranging from 20 to 80 L/min at FiO<sub>2</sub> ranging from 21% to 90%. Testing was performed on heated and humidified air.

## Results

There were significant differences between the actual flow and the set flow in all five devices ( $p < 0.001$ ). The AIRVO 2 had the lowest differences followed by the bellavista 1000. TNI softFlow compensated the best when one of the two prongs of the nasal cannula was obstructed. There were significant differences in actual FiO<sub>2</sub> compared to set FiO<sub>2</sub> in all devices. The AIRVO 2 showed the least difference, followed by the bellavista™ 1000. In devices other than the bellavista™ 1000, increasing the set FiO<sub>2</sub> resulted in larger differences between set and measured FiO<sub>2</sub>.



**Figure 1** Flow deviation at different set-flow levels in five devices

## Why use the bellavista™ for HFNC therapy?

HFNC oxygen therapy offers significant physiologic benefits for patients requiring respiratory support without ventilation. HFNC-therapy has multiple physiological advantages, including generating positive end-expiratory pressure (PEEP), washing out dead space, decreasing inspiratory effort, improving alveolar recruitment, and increasing lung volume and compliance. Using a HFNC system allows the consistent delivery of warmed, humidified oxygen therapy to the patient who is struggling to breathe. Wall oxygen is dry and cold and far less consistent in flow and  $\text{FiO}_2$  delivered. The systems tested in this study offer distinct therapeutic benefits compared to conventional oxygen therapy.

Not every patient succeeds with HFNC, with some requiring invasive or noninvasive ventilation. Being able to provide HFNC therapy using a ventilator is a potential game-changer for physicians. If a patient requires additional support such as invasive or noninvasive ventilation, the life-saving machine is already at the patient's bedside. The bellavista 1000

ventilator can provide high quality and predictable HFNC therapy without requiring the capital outlay of a stand-alone machine.



Figure 2 The bellavista™ 1000 ventilator

## Take home message

According to this bench study, the bellavista™ 1000 ventilator with the HFNC module delivers high flow nasal cannula oxygen therapy at a level comparable to the best dedicated free-standing products on the market. Flow rates and  $\text{FiO}_2$  are consistent and predictable even in cases of accidental compression or blocking of the cannulas. The bellavista has the additional advantage of being able to be used for higher levels of ventilatory support should the patient require it.

### GLOBAL HEADQUARTERS

Vyaire Medical, Inc. | 26125 North Riverwoods Blvd | Mettawa, IL 60045, USA

#### For global distribution.

© 2021 Vyaire. Vyaire, the Vyaire logo and all other trademarks or registered trademarks are property of Vyaire Medical, Inc., or one of its affiliates. VYR-GBL-2100179