

Subject: SentrySuite™ – SQL Injection and Authenticated Remote Code Execution vulnerabilities

Date: 2022-12-20 (last update)

Originator: Timo Kosig

Background

Vyairé proactively communicates with our customers to inform them about cybersecurity issues to help healthcare delivery organizations identify and remediate potential cybersecurity risks.

This bulletin contains product security information and recommendations related to two cybersecurity vulnerabilities that Vyairé identified within its SentrySuite™ software.

Response

Affected Vyairé Products

Vyairé has confirmed that the following versions of its SentrySuite™ software are affected by the vulnerabilities found:

- SentrySuite™ 3.20 (up to and including Maintenance Package 5)

Vyairé is currently reviewing the status of the following SentrySuite versions which are still under active support by Vyairé:

- SentrySuite™ 3.10
- SentrySuite™ 3.0

This bulletin will be updated with further information as it becomes available.

Product Security Bulletin

Important Information – Please Read and Keep



Affected Product Components

Vyairé SentrySuite™ software offers two different deployment options, SentrySuite™ can either be deployed as a client-server application with a dedicated application backend server, or as a stand-alone application.

Please see the following table to identify which product components are affected depending on the chosen deployment scenario.

Deployment mode	Vulnerability	
	SQL Injection	Authenticated Remote Code Execution
Stand-alone	Only if connected to the network	Only if connected to the network
Client-server	Yes, only application backend server is affected, clients are not affected	Yes, only application backend server is affected, clients are not affected

Vulnerability Details

SQL Injection

The SentrySuite™ application backend server and SentrySuite™ installed as a stand-alone system are vulnerable to a SQL Injection attack. Vyairé Product Security has analyzed the vulnerability and has assigned the following CVSS score to this vulnerability:

CVSS: 8.3 (High) CVSS:3.1/AV:N/AC:L/PR:L/UI:N/S:U/C:H/I:H/A:H

Authenticated Remote Code Execution

The SentrySuite™ application backend server and SentrySuite™ installed as a stand-alone system are vulnerable to an Authenticated Remote Code Execution attack. Vyairé Product Security has analyzed the vulnerability and has assigned the following CVSS score to this vulnerability:

CVSS: 8.2 (High) CVSS:3.1/AV:N/AC:L/PR:L/UI:N/S:U/C:H/I:H/A:H

Mitigations & Compensating Controls

Vyair recommends implementing the following mitigations and compensating controls to reduce risk associated with these vulnerabilities:

For customers running SentrySuite[™] 3.20 (up to and including Maintenance Package 5)

- Install SentrySuite[™] 3.20 Maintenance Package 6 which fully remediates both vulnerabilities. Please contact your Vyair field service representative or get in touch via the [Vyair website](#) to schedule installation.

For customers running SentrySuite[™] 3.10 or 3.0

- For client-server deployments
 - Please follow the generic controls described in the next section.
- For network-connected stand-alone deployments:
 - Ensure that the host-based firewall is enabled on the computer system and is configured correctly for your usage. SentrySuite[™] does not require any open incoming ports in stand-alone mode in a default configuration. Usage of electronic interfaces to hospital information systems, electronic medical record systems or other systems must be evaluated separately.

Product Security Bulletin

Important Information – Please Read and Keep



Generic controls that should be considered in all scenarios

The following generic controls are recommended by the document *Vyairé Product Security White Paper - SentrySuite 3.20 V1* which is available to all customers upon request:

- Ensure your data has been backed up and stored according to your individual process and that your disaster recovery procedures are in place.
- Update your anti-virus and malware protection, where available.

For product or site-specific concerns, contact your Vyairé service representative.

For more information on Vyairé's proactive approach to product security and vulnerability management, contact us at productsecurity@vyaire.com or visit www.vyaire.com/product-security

A handwritten signature in black ink, appearing to read "Timo Kosig".

Timo Kosig
timo.kosig@vyaire.com
Product Security Manager
Vyairé Medical GmbH

A handwritten signature in black ink, appearing to read "Holger Schmitt".

Holger Schmitt
holger.schmitt@vyaire.com
Senior Product Manager Software and IT integration
Vyairé Medical GmbH

GLOBAL HEADQUARTERS

Vyairé Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA

26125 N. Riverwoods Blvd., Mettawa, IL 60045, USA

vyaire.com For global distribution.

Trademarks are the property of their respective owners.

© 2022 Vyairé. Vyairé and the Vyairé logo and all other trademarks or registered trademarks are property of Vyairé Medical, Inc., or one of its affiliates. Please read the complete Instructions for Use that come with the devices or follow the instructions on the product labelling. | VYR-GBL-2200190 (1.0)