

vyaire™
M E D I C A L



Vyaire AirLife™
Noninvasive ventilation masks



Noninvasive ventilation masks

Our noninvasive ventilation (NIV) full face and nasal masks offer clinicians a reliable option, while providing patients comfort and security. The adjustable four-point headgear and flexible swivel connector enable ease in fit and adjustment.

Masks feature a lightweight silicone design that contours to the patient's face for a comfortable fit, while the quick release clip eliminates the need to refit after removal. The odorless, easy-to-clean masks include both vented and non-vented options and are available in three sizes to fit a wide range of patients.

NIV masks (10/cs.)

Cat no.	Product description
NIV040S	Full face mask, vented with anti-asphyxia valve (AAV), small
NIV040M	Full face mask, vented with AAV, medium
NIV040L	Full face mask, vented with AAV, large
NIV043S	Full face mask, non-vented with AAV, small
NIV043M	Full face mask, non-vented with AAV, medium
NIV043L	Full face mask, non-vented with AAV, large
NIV041S	Full face mask, non-vented, small
NIV041M	Full face mask, non-vented, medium
NIV041L	Full face mask, non-vented, large
NIV042S	Nasal mask, vented, small
NIV042M	Nasal mask, vented, medium
NIV042L	Nasal mask, vented, large
NIV044S	Nasal mask, non-vented with AAV, small
NIV044M	Nasal mask, non-vented with AAV, medium
NIV044L	Nasal mask, non-vented with AAV, large

Key features

- Flexible swivel connector enables ease in positioning of the circuit
- Comfortable four-point headgear with crown strap for better stability and sealing
- Adjustable forehead rest for comfort and stability
- Soft seal that conforms to the face, reducing need to overtighten
- Quick release strap
- Made without DEHP

NIV nasal masks



Quick release strap



Adjustable forehead rest





For more information, please contact your local Vyairē representative.

GLOBAL HEADQUARTERS

Vyairē Medical, Inc.
26125 N. Riverwoods Blvd.
Mettawa, IL 60045
USA



EMERGO EUROPE
Prinsessegracht 20
2514 AP The Hague
The Netherlands



Vyairē Medical, Inc.
26125 North Riverwoods Blvd
Mettawa, IL 60045
USA



AUSTRALIAN SPONSOR

Vyairē Medical Pty Ltd
Level 5, 7 Eden Park Drive
Macquarie Park, NSW, 2113
Australia

vyairē.com

For Europe, Australia, Asia, Middle East and Africa distribution only.

© 2019 Vyairē. Vyairē, the Vyairē logo and AirLife are trademarks or registered trademarks of Vyairē Medical, Inc., or one of its affiliates. Medical devices class IIa according to Medical Devices Directive 93/42/EC. Please read the complete Instructions For Use that come with the devices or follow the instructions on the product labelling. VYR-INT-1900086